What's NEW in Land 1985.

Automatic 360° Body Stitching for XT and X3 captures Simply take a front and back view of the patient and VECTRA software will automatically combine them into a single three-dimensional wrap-around image.

Dynamic photo mode

adds a third shading mode, similar to Gray mode. The image is skin colored, but with reduced detail.



anatomical features under raked lighting and accentuates cleavage.

This helps

shape and

to show

NEW in Breast Sculptor®

- Improved Mastopexy simulations More breast shapes are now accommodated.
- Breast implant libraries are updated with the latest changes made by the implant manufacturers.

Animation tool now enables you to incorporate multiple images into a single animation, such as fading between pre-op and simulated images, or pre-op and post-op images.

Combine pre-op and post-op images to create animations for your website



Windows 10 compatible for all consultation computers and VECTRA H1 capture computers.

IMAGING EXCELLENCE FROM



www.canfieldsci.com / info@canfieldsci.com / phone +1.973.276.0336 / (USA) 800.815.4330 3D SOLUTIONS · FACIAL IMAGING & ANALYSIS · IMAGING SOFTWARE · PHOTOGRAPHY · RESEARCH SYSTEMS & SERVICES · TRAINING